

Aware (QS-110P)

The Aware™ Alarm Management Server Extension Captures alarm information by type and device. It also captures the reason for alarm, i.e., lower threshold or higher threshold violated. This information can be utilized by hospitals to tweak alarm settings to avoid nuisance alarms. Additionally, Aware™ goes beyond simple aggregation and provides meaningful and actionable patient centric reports. Finally, the data is stored in a sequel database for use during report generation.

Hardware/Software requirements:

1. 64-bit Windows 7 or 64-bit Windows 10
2. Microsoft .NET Framework 4.5.2 (automatically installed)
3. Microsoft Visual C++ 2015 Redistributable (automatically installed)
4. 3 GB RAM
5. 1 GB of available hard disk space

Aware plug-in:

- It is running inside of Enterprise Gateway (US Version).

Length of storage:

- 62 days
- 200 million records of alarms

* Alarms older than 62 days will be deleted automatically.

