

## Digital Health Platform (1010)

The NKDHS Digital Health Platform stores high fidelity data (vital signs, alarms, and waveforms) from Nihon Kohden patient monitors (and compatible third-party devices). The system includes a web-based application, where data can be retroactively viewed and extracted. The Digital Health Platform also support storage patient demographic data from an ADT feed.

### Cloud infrastructure supported:

1. Amazon Web Services
2. Dell Virtustream

\* We are cloud provider agnostic.

### Streaming is composed of 4 servers:

	<b>CPU Cores</b>	<b>Memory</b>	<b>Storage</b>	<b>Network</b>	<b>OS</b>
A	16	32	200 GB	1 Gigabit	Windows Server 2016 Standard
B	16	32	500 GB	1 Gigabit	Ubuntu Desktop/Server 18.04.4 LTS
C	32	64	10 TB	1 Gigabit	Windows Server 2016 Standard
D	16	32	200 GB	1 Gigabit	Windows Server 2016 Standard

\* 10 TB- main storage

### Product image:

