

[NK-HIQ-Enterprise Gateway International \(QP-988P-INT\)](#)**SUPPORT:****Application supported:**

1. ECG Server Extension (QP-992P)
2. HL7 Server Extension (CGS-9012)
3. Wireless Gateway (QP-990P)
4. ViTrac (QP-989P)
5. Pager Gateway (CGS-9010)

Max number for EGI are:

1. Clients: 100
2. Wi-Fi Devices: 340
3. Total Bed: 500

SYSTEM REQUIREMENTS:**Software Requirements:**

1. Window Server 2019 Standard
2. Microsoft SQL Server 2019 Express
3. Strawberry Perl

Hardware requirement:

1. 3.3 GHz 4 Core/4 Thread Server Class Processor
* Example: 3.3 GHz Intel Xeon E-2124 Processor (or better)
2. 8 GB RAM
3. 1 TB 7200RPM Hard Drive
4. 2 Gigabit Network Interface Controller

*Recommended hardware Dell R240.

Product Image: