

NK-HiQ Enterprise Gateway (QP-988P-A)**Requirements:****Hardware requirements:**

1. Dell R640
2. 2.3 GHz Intel Xeon Gold 5118 Processor
3. 64 GB RAM
4. 2 1.2TB 10K RPM SAS Hard Drive (RAID1)
5. 3 Gigabit Network Interface Controllers
6. 1 Internal USB Port

Software requirements:

1. Windows Server 2019 Standard
2. Microsoft SQL Server 2019 Standard
3. Strawberry Perl

Support:**Applications supported:**

1. Wireless Gateway
2. ViTrac
3. RNS
4. HL7
5. ECG
6. Pager
7. Waveform
8. Aware

Max beds supported:

- 2,500

Max clients supported:

- 300

Max Wi-Fi Devices supported:

- 1,500

Product Image: