

NKAnywhere (QP-999P)

NKAnywhere™ allows for software downloads and offers the ability to remotely set settings for Nihon Kohden devices. Additionally, NKAnywhere can execute multiple installations on multiple devices concurrently, then it monitors and reports the installation progress in the user interface screens.

Requirements:

* NKAnywhere requires sufficient processing power and memory to monitor and install software to multiple Agent-enabled devices; the application itself does not consume a lot of processing power.

NKAnywhere admin:

- Windows 10
- 2GB of memory or higher
- 2GB of available hard disk space or more

Agent:

* NKAnywhere Windows Agent application runs on devices with

Windows Operating system:


- Windows XP, 7 and 10
- 2GB of memory or higher
- 1GB of available hard disk space

Support:**Devices supported:**

See table on page 2.

Concurrent connections supported:

- Admin application allows installing software on 20 wired or 5 wireless devices concurrently.



System	Agent Type	Models	Supported devices for firmware update	Supported devices for setting update
BSM	Device	BSM-1700 Series	Ver. 01-15 or later.	Ver. 02-10 or later.
		BSM-3000 Series	Ver. 07-02 or later.	Ver. 07-41 or later.
		BSM-6000 Series	Ver. 07-02 or later.	Ver. 07-41 or later.
CSM	Device	CSM-1901	Ver. 01-21 or later.	Ver. 01-27 or later.
ORG	Device	ORG-9100 series	Ver. 04-01 or later.	N/A.
GZ	Device	GZ-120P	Ver. 02-01 or later.	Ver. 02-01 or later.
		GZ-130P		
		GZ-140P		
SVM	Device	SVM-7200 Series	Ver. 02-02 or later	Ver. 02-02 or later
CNS	Windows	CNS-6201	Ver. 05-10 or later.	Ver. 05-10 or later.
		CNS-6801	Ver. 02-10 or later.	Ver. 02-10 or later.
CRC	Windows	RNS-6803	Ver. 02-10 or later	Ver. 02-10 or later
		RNS-6804		