

Ventilator Gateway (CGS-9013)

The Ventilator Gateway is a software application that relays data from Nihon Kohden OrangeMed ventilators to a hospital information system (HIS). The Gateway can be installed and configured to expand clinician access to patient data.

Support:

- 150 Ventilators per installation

Compatibility:

- NKV-550 Ventilator
- NKV-330 Ventilator

Hardware/ Software Requirements:

- I. Running on a 64-bit version of Windows
- II. Proper licensing for the Ventilator interface is required
- III. Communicates over a wireless network

Product image: