

ViTrac (QP-989P)

ViTrac™ is designed to securely provide Physicians, Nurses, and other medical personnel with a remote display of medical device data from the Nihon Kohden Central Nursing Station (CNS), that includes a wide range of vital signs on multiple patients. The data can be viewed on Mobile devices within the hospital network or remotely, via a VPN connection, on-demand and in near real-time. ViTrac is available on iOS, Android, and Windows devices.

**Requirements:****Hardware requirements (WINDOWS):**

- Recommended screen resolution: 1920 x 1200  
(Minimum screen resolution: 1920 x 1080)
- Aspect ratio 16:9
- 8<sup>th</sup> Gen Intel Core, 4 CPU Cores
- 16 GB RAM
- 256 GB Hard Disk

**Software requirements:**

- Windows 10

**Support:****Number of beds supported:**

- Tablet= 12
- Phone= 8

**IOS Version supported:**

- Devices running iOS 11 or newer

**Android Version supported:**

- Devices running Android OS 9, 10, or 11

**EG applications needed:**

- Compatible with both International EG and US EG versions.

**Product image:**